



# Update to EN 61010-1 Amendment Statement of Operating Conditions

Catalog #	Description
1652960	Mini-PROTEAN 2-D Cell
1652961	Mini-PROTEAN Tube Cell
1652965	Mini-PROTEAN Tube Cell Module
1703930	Mini Trans-Blot Electrophoretic Transfer Cell
1703935	Mini Trans-Blot Module
1703989	Mini Trans-Blot Cell and PowerPac Basic Power Supply
1703836	Mini Trans-Blot Cell and PowerPac HC Power Supply
1652976	Model 422 Electro-Eluter
1704070	Criterion Blotter with Plate Electrodes
1704071	Criterion Blotter with Wire Electrodes
1656001	Criterion Cell

---

For research purposes only.

---

## Specifications

See the following user manual updates according to the guidelines in IEC/EN/UL 61010-1/A1, section 5.4.2.

## Environmental Conditions

Operating temperature	0–50°C, indoor use
Humidity	10–95%, in the absence of condensation
Altitude	Up to 2,000 m
Pollution degree	2

## Input Power

Supply voltage fluctuations	Not applicable
Overvoltage category	Not applicable

BIO-RAD, MINI-PROTEAN, and TRANS-BLOT are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.