



Update to EN 61010-1 Amendment Statement of Operating Conditions

Catalog #	Description
1703940	Trans-Blot SD Semi-Dry_Electrophoretic Transfer Cell

For research purposes only.

Specifications

See the following user manual updates according to the guidelines in IEC/EN/UL 61010-1/A1, section 5.4.2.

Environmental Conditions

Operating temperature	0–50°C, indoor use
Humidity	10–95%, in the absence of condensation
Altitude	Up to 2,000 m
Pollution degree	2

Input Power

Supply voltage fluctuations	Not applicable Maximum operating voltage: DC 25 V Maximum power: 50 W
Overvoltage category	Not applicable

BIO-RAD and TRANS-BLOT are trademarks of Bio-Rad Laboratories, Inc. in certain jurisdictions. All trademarks used herein are the property of their respective owner.